



# REBOUND™

Introducing Rebound **Full Thickness** Placental-Derived Allograft **Matrix.**

**Rebound™ consists of all 3 layers of the placental membrane and is intended for use as a barrier or wound cover.**

A full thickness placental derived allograft matrix to compliment the Celularity product portfolio.

- ◆ Dehydrated and decellularized
- ◆ Terminally sterilized
- ◆ Flexible to use in a variety of wounds
- ◆ Store at ambient room temperature 15-30°C
- ◆ 5-year shelf life

## CONTRAINDICATIONS

Contraindications for the use of this allograft include, but are not limited to:

- ① The presence of infection at the surgical site and/or distant foci of infections that may spread to the surgical site;
- ② Uncooperative patient or a patient with neurologic disorders who is incapable of following direction, including weight control and activity levels;
- ③ Pregnancy.

## POSSIBLE ADVERSE EFFECTS

Any adverse outcomes potentially attributable to this graft must be reported immediately to Sequence at [Regulatory@sequencelifesci.com](mailto:Regulatory@sequencelifesci.com)

Available in 6 different **sizes** to meet a variety of wounds



Product Code:  
REB27248



Product Code:  
REB27246



Product Code:  
REB27244



Product Code:  
REB27224



Product Code:  
REB27223



Product Code:  
REB27222

**REBOUND**<sup>TM</sup>

**For more information or to place an order  
please contact Celularity at  
1-844-963-2273**